SIEMENS

Siemens Healthcare Diagnostics Reagent Strips for Urinalysis

SUMMARY AND EXPLANATION: Siemens Healthcare Diagnostics Reagent Strips for Urinalysis include test pads for protein, blood, leucocytes, nitrite, glucose, ketone (acetoacetic acid), pH, specific gravity, bilirubin, and urobilinogen. Refer to the carton or bottle label for the tests included on the product you are using. The strips are for professional, *in vitro* diagnostic use (IVD) only. Read the insert carefully before using the product ().


Siemens Reagent Strips are ready to use upon removal from the bottle. The strips may be read visually. They can also be read instrumentally, using the CLINITEK® family of Urine Chemistry Analysers and the appropriate software; contact your product representative for further information.

CAUTION: Ensure that work areas and specimen containers are always free of detergents and other contaminating substances. Some substances can interfere with patient results.

SPECIMEN COLLECTION AND PREPARATION: Collect freshly-voided urine in a clean, dry container. Mix the sample before testing and test it within two hours after voiding, sooner if testing for bilirubin or urobilinogen. Contamination of the urine specimen with skin cleansers containing chlorhexidine may affect protein (and to a lesser extent specific gravity and bilirubin) test results. Work areas and specimen containers should always be free of detergents and other contaminating substances. If unable to test within the recommended time, refrigerate the specimen immediately and let it return to room temperature before testing.

VISUAL PROCEDURE:

1. Dip all the test pads of the strip into the urine and remove immediately.
2. Drag the edge of the strip against the container rim to remove excess urine.
3. Compare each test pad to the colour blocks on the bottle label. Read each pad at the time shown on the label, starting with the shortest time. Colour changes that occur after 2 minutes are of no diagnostic value. Discard the used reagent strip according to your laboratory procedures.

INSTRUMENTAL PROCEDURE: Carefully follow the directions given in the appropriate instrument operating manual.

QUALITY CONTROL: Test known negative and positive specimens or controls whenever a new bottle is first opened. Water should NOT be used as a negative control. Each laboratory should establish its own goals for adequate standards of performance. CHEK-STIX® Positive and Negative Control Strips provide a convenient basis for a quality control programme.

STORAGE AND HANDLING: Store at temperatures between 15°–30°C (59°–86°F). Do not use the strips after their expiration date (). Do not store the bottle in direct sunlight and do not remove the desiccant from the bottle. PROTECTION AGAINST EXPOSURE TO LIGHT, HEAT AND AMBIENT MOISTURE IS MANDATORY TO GUARD AGAINST ALTERED REAGENT REACTIVITY. Do not remove the strip from the bottle until immediately before it is to be used for testing. Replace the cap immediately and tightly after removing the reagent strip. Do not touch the test areas of the strip. Discolouration or darkening of the test pads may indicate deterioration. If this is evident, or if test results are questionable or inconsistent with expected findings, confirm that the product is within its expiration date and is reacting properly using known negative and positive control materials.

LIMITATIONS OF PROCEDURE: As with all laboratory tests, definitive diagnostic or therapeutic decisions should not be based on any single result or method. Substances that cause abnormal urine colour may affect the readability of test pads on urinalysis reagent strips. These substances include visible levels of blood or bilirubin and drugs containing dyes, nitrofurantoin, or riboflavin. Levels of ascorbic acid normally found in urine do not interfere with these tests.

TEST INFORMATION: **PROTEIN:** Less than 0.15 g (150 mg) of total protein is normally excreted per day (24 hour period). Clinical proteinuria is indicated at greater than 0.5 g (500 mg) of protein per day

(strip result of ≥ 0.3 g/L or 30 mg/dL). Clinical judgement is needed to evaluate the significance of Trace results. The protein test is less sensitive to mucoproteins and globulins, which are generally detected at levels of 0.6 g/L (60 mg/dL) or higher; a negative result does not rule out the presence of these other proteins.

BLOOD: Normally, no haemoglobin is detectable in urine (< 100 μ g/L or 0.010 mg/dL; 3 RBC/ μ L). The significance of the Trace reaction may vary among patients, and clinical judgment is required for assessment in an individual case. Blood is often, but not always, found in the urine of menstruating females. The test is equally sensitive to myoglobin as to haemoglobin. A haemoglobin concentration of 150–620 μ g/L (0.015–0.062 mg/dL) is approximately equivalent to 5–20 intact red blood cells per microliter. Captopril and other compounds that contain sulphydryl groups may reduce the sensitivity. Certain oxidizing contaminants, such as hypochlorite, may produce false positive results. Microbial peroxidase associated with urinary tract infection may cause a false positive reaction.

LEUCOCYTES: Normal urine specimens generally yield negative results. A strip result of Small or greater is a useful indicator of infection. Trace results may be of questionable clinical significance; however, Trace results observed repeatedly may be clinically significant. Elevated glucose concentrations (≥ 160 mmol/L or 3 g/dL) may cause decreased test results. The presence of cephalixin, cephalothin, or high concentrations of oxalic acid may also cause decreased test results. Tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. Positive results may occasionally be due to contamination of the specimen by vaginal discharge.

NITRITE: Normally, no nitrite is detectable in urine. This test depends upon the conversion of nitrate (derived from the diet) to nitrite by the action of Gram negative bacteria in the urine. Many enteric gram-negative organisms give positive results when their number is greater than 10^5 /mL (16.2 μ mol/L or 0.075 mg/dL nitrite ion or greater). The test is specific for nitrite and will not react with any other substance normally excreted in urine. Pink spots or pink edges should not be interpreted as a positive result. A negative result does not rule out significant bacteriuria. False negative results may occur with shortened bladder incubation of the urine (< 4 hours), absence of dietary nitrate, or the presence of nonreductive pathological microbes.

GLUCOSE: Small amounts of glucose (< 1.67 mmol/L or 30 mg/dL) are normally excreted by the kidney. These amounts are usually below the sensitivity level of this test but on occasion may produce a result between Negative and 5.5 mmol/L (100 mg/dL) that is interpreted as a positive result. The test is specific for glucose; no substance excreted in urine other than glucose is known to give a positive result. High ketone levels (4 mmol/L or 40 mg/dL) may cause false negatives for specimens containing small amounts of glucose (4–7 mmol/L or 75–125 mg/dL).

KETONE: Normally, no ketone is detectable in urine. The test reacts with acetoacetic acid in urine. It does not react with acetone or β -hydroxybutyric acid. False Trace results may occur with highly pigmented urine specimens or those containing large amounts of levodopa metabolites. Compounds that contain sulphydryl groups, such as mesna (2-mercaptoethane sulfonic acid) and captopril, may cause false positive results or an atypical colour reaction.

pH: The normal pH of urine can range from 4.6 to 8.0. The pH test area measures pH values from 5–8.5 visually and 5–9 instrumentally, generally to within one unit of the expected result. Bacterial growth by certain organisms in a specimen may cause a marked alkaline shift (pH > 8.0), usually because of urea conversion to ammonia.

SPECIFIC GRAVITY: The normal SG of urine ranges from 1.001 to 1.035. If the SG of a random urine specimen is ≥ 1.023 , the concentrating ability of the kidneys can be considered normal. This test permits determination of urine specific gravity between 1.000 and 1.030. In general, it correlates within 0.005 with values obtained with the refractive index method. For increased accuracy, 0.005 may be added to visual readings from urines with pH ≥ 6.5 . Strips read instrumentally are automatically adjusted for pH by the instrument. The Siemens SG test is not affected by the presence of radiopaque dyes as are the refractive index, urinometer, and osmolality methods. Highly buffered alkaline urines may cause low readings, while the presence of moderate quantities of protein (1–7.5 g/L or 100–750 mg/dL) may cause elevated readings.

BILIRUBIN: Normally, no bilirubin is detectable by even the most sensitive methods. Even trace amounts of bilirubin are sufficiently abnormal to require further investigation. Indican (indoxyl sulfate) can produce a yellow-orange to red colour response that may interfere with the interpretation of a negative or positive reading. Metabolites of etodolac may cause false positive or atypical results. Atypical colours may indicate bile pigment abnormalities and the urine specimen should be tested further.

UROBILINOGEN: Urobilinogen is normally present in urine at concentrations up to 16 μ mol/L (1.0 mg/dL). A result of 33 μ mol/L (2.0 mg/dL) represents the transition from normal to abnormal, and the patient and/or urine specimen should be evaluated further. This test area will detect urobilinogen in concentrations as low as 3.2 μ mol/L (0.2 mg/dL or 0.2 EU/dL) in urine. The absence of urobilinogen in the specimen cannot be determined. The test pad may react with interfering substances known to react with Ehrlich's reagent, such as β -aminosalicylic acid and sulfonamides. Atypical colour reactions may be obtained in the presence of high concentrations of β -aminobenzoic acid. False negative results may be obtained if formalin is present. Strip reactivity increases with temperature; the optimum temperature is 22°–26°C. The test is not a reliable method for the detection of porphobilinogen.

SPECIFIC PERFORMANCE CHARACTERISTICS: Performance characteristics are based on clinical and analytical studies and depend upon several factors: the variability of

colour perception; the presence or absence of inhibitory and matrix factors typically found in urine; and the laboratory conditions in which the product is used (e.g., lighting, temperature, and humidity). Each colour block or instrumental result represents a range of values. Because of specimen and reading variability, specimens with analyte concentrations that fall between nominal levels may give results at either level. Results will usually be within one level of the true concentration. Exact agreement between visual results and instrumental results might not be found because of the inherent differences between the perception of the human eye and the optical systems of the instruments. The following list shows the generally detectable levels of the analytes in contrived urines; however, because of the inherent variability of clinical urines, lesser concentrations may be detected under certain conditions.

Test Pad and Sensitivity:

Protein: 0.15–0.3 g/L (15–30 mg/dL) albumin

Blood: 150–620 μ g/L (0.015–0.062 mg/dL) haemoglobin

Leucocytes: 5–15 cells/hpf in clinical urine

Nitrite: 13–22 μ mol/L (0.06–0.1 mg/dL) nitrite ion

Glucose: 4–7 mmol/L (75–125 mg/dL) acetoacetic acid

Ketone: 0.5–1.0 mmol/L (5–10 mg/dL) acetoacetic acid

Bilirubin: 7–14 μ mol/L (0.4–0.8 mg/dL) bilirubin

CHEMICAL PRINCIPLES OF PROCEDURES AND INGREDIENTS: (based on dry weight at time of impregnation)

Protein: This test is based on the protein-error-of-indicators principle. **Ingredients:** 0.3% w/w tetrabromophenol blue; 97.3% w/w buffer; 2.4% w/w nonreactive ingredients.

Blood: This test is based on the peroxidase-like activity of hemoglobin, which catalyzes the reaction of diisopropylbenzene dihydroperoxide and 3,3',5,5'-tetramethylbenzidine. **Ingredients:** 6.8% w/w diisopropylbenzene dihydroperoxide; 4.0% w/w 3,3',5,5'-tetramethylbenzidine; 48.0% w/w buffer; 41.2% w/w nonreactive ingredients.

Leucocytes: Granulocytic leucocytes contain esterases that catalyze the hydrolysis of the derivatized pyrrole amino acid ester to liberate 3-hydroxy-5-phenyl pyrrole. This pyrrole then reacts with a diazonium salt. **Ingredients:** 0.4% w/w derivatized pyrrole amino acid ester; 0.2% w/w diazonium salt; 40.9% w/w buffer; 58.5% w/w nonreactive ingredients.

Nitrite: At the acid pH of the test pad, nitrite in the urine reacts with ρ -arsanilic acid to form a diazonium compound. This diazonium compound in turn couples with 1,2,3,4-tetrahydrobenzo(h)quinolin-3-ol. **Ingredients:** 1.4% w/w ρ -arsanilic acid; 1.3% w/w 1,2,3,4-tetrahydrobenzo(h)quinolin-3-ol; 10.8% w/w buffer; 86.5% w/w nonreactive ingredients.

Glucose: This test is based on a double sequential enzyme reaction. Glucose oxidase catalyzes the formation of gluconic acid and hydrogen peroxide from the oxidation of glucose. Peroxidase then catalyzes the reaction of hydrogen peroxide with a potassium iodide chromogen to oxidize the chromogen. **Ingredients:** 2.2% w/w glucose oxidase (microbial, 1.3 IU); 1.0% w/w peroxidase (horse radish, 3300 IU); 8.1% w/w potassium iodide; 69.8% w/w buffer; 18.9% w/w nonreactive ingredients.

Ketone: This test is based on the development of colours when acetoacetic acid reacts with nitroprusside. **Ingredients:** 7.1% w/w sodium nitroprusside; 92.9% w/w buffer.

pH: This test is based on a double indicator principle that gives a broad range of colours covering the entire urinary pH range. **Ingredients:** 0.2% w/w methyl red; 2.8% w/w bromthymol blue; 97.0% w/w nonreactive ingredients.

Specific Gravity: This test is based on the apparent pKa change of certain pretreated poly-electrolytes in relation to ionic concentration. **Ingredients:** 2.8% w/w bromthymol blue; 68.8% w/w poly (methyl vinyl ether/maleic anhydride); 28.4% w/w sodium hydroxide.

Bilirubin: This test is based on the coupling of bilirubin with diazotized dichloroaniline in a strongly acid medium. **Ingredients:** 0.4% w/w 2,4-dichloroaniline diazonium salt; 37.3% w/w buffer; 62.3% w/w nonreactive ingredients.

Urobilinogen: This test is based on the Ehrlich reaction in which ρ -diethylamino-benzaldehyde in conjunction with a colour enhancer reacts with urobilinogen in a strongly acid medium. **Ingredients:** 0.2% w/w ρ -diethylaminobenzaldehyde; 99.8% w/w non-reactive ingredients.

US PATENT NUMBERS: Refer to the carton of the product you are using for the US patent numbers.

TRADEMARKS: MULTISTIX®, CLINITEK®, and CHEK-STIX® are trademarks of Siemens Healthcare Diagnostics.

PRODUCT NOS.: 2289, 2300, 2304, 2308, 2740, 2741, 2743, 2810, 2815, 2820, 2857, 2877.

For more information, contact your Siemens representative or Customer Service.

Siemens Healthcare Diagnostics Inc.
Sir William Siemens Square
Frimley, Camberley
Surrey, GU16 8QD UK
Direct Telephone: +44 (0) 845 600 1955
Direct Fax: +44 (0) 1276 696 680

SIEMENS

Siemens Healthcare Diagnostics Inc.
Tarrytown, NY 10591-5097 USA

Siemens Healthcare Diagnostics Ltd.
Sir William Siemens Sq.
Frimley, Camberley, UK GU16 8QD

www.siemens.com/diagnostics

© 2008 Siemens Healthcare Diagnostics. All rights reserved.

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