

schülke -+



Active oxygen-based concentrate in powder form for disinfecting and cleaning medical devices and surfaces.

perform®



Our Plus

- very good material compatibility
- broad spectrum of activity (limited spectrum of virucidal activity, tuberculocidal, sporicidal (incl. *C. difficile*))
- active oxygen based

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal EN16615 - clean conditions	0.5 % (5 g/l)	60 min.
bactericidal EN13727 - clean conditions	0.5 % (5 g/l)	5 min.
bactericidal EN16615 - clean conditions	1 % (10 g/l)	15 min.
bactericidal EN13727 - clean conditions	1 % (10 g/l)	1 min.
tuberculocidal EN14348 - clean conditions	2 % (20 g/l)	60 min.
yeastcidal EN16615 - clean conditions	0.5 % (5 g/l)	60 min.
yeastcidal EN13624 - clean conditions	0.5 % (5 g/l)	60 min.
yeastcidal EN16615 - clean conditions	1 % (10 g/l)	15 min.
yeastcidal EN13624 - clean conditions	1 % (10 g/l)	15 min.
fungicidal EN13624 - dirty conditions	2 % (20 g/l)	60 min.
virucidal against enveloped viruses EN14476	0.5 % (5 g/l)	5 min.

Efficacy	Concentration	Contact time
limited spectrum virucidal activity EN14476 - clean conditions	0.5 % (5 g/l)	5 min.
virucidal EN14476 - clean conditions	1 % (10 g/l)	5 min.
Clostridium difficile EN17126 - clean conditions	1 % (10 g/l)	2 h
Clostridium difficile EN13697 - clean conditions	1 % (10 g/l)	2 h
Clostridium difficile EN17126 - clean conditions	2 % (20 g/l)	60 min.

Application area	Concentration	Contact time
RKI-listed according to German law (§18 IfSG) - scope AB	3 % (30 g/l)	4 h

Certificates

- RKI listed
- IHO listed
- VAH certificate



Product data

Composition:

100 g powder contains the following active substances: 45 g

Pentapotassium-bis-(peroxymonosulphate) bis(sulphate)

Labelling according to Regulation (EC) No. 648/2004:

5 - 15 % anionic surfactants; < 5 % non-ionic surfactants; < 5 % soap and phosphonates, perfumes

Chemical-physical data

Color	white
Flash point	Not applicable
Form	granular
pH	ca. 4 / 5 g/l / 20 °C / in water
Viscosity, dynamic	Not applicable

Special advice

Use disinfectants safely. Always read the label and product information before use.

The working solution may have a bleaching effect on certain very sensitive textile dyestuffs. Wearing protective gloves is advised. Store the container tightly closed and dry. For further information on safe handling and storage conditions, please refer to the safety data sheet. The shelf-life of the opened container can be found on the packaging. Ensure adequate ventilation at the site of application. Do not store bags or boxes above room temperature.

Application areas

perform® is suitable for disinfecting surfaces and medical devices. As a very well established product with a broad spectrum of activity and high material compatibility, it is used in all hospital areas. perform is also suitable for disinfecting respiratory and inhalation masks as well as training dummies.

Instructions for use

To prepare the working solution, fill the bucket with water, add the product and dissolve. When using the box, the corresponding dosing is carried out according to the dosing table or the box label:

- 1 bag (~40 g) to 8 l of water = 0.5% solution

- 2 bags (~80 g) to 8 l of water = 1% solution

The prepared working solutions have a lifetime of max. 30 hours (at 20°C) and 4 hours (at 40°C) for guaranteed stability of the active ingredient. The working solution should be replaced earlier in case of heavy contamination. Apply the working solution in the desired concentration to the floor using the 2-step wiping procedure. Wet wipe medical devices and other surfaces. To avoid formation of residues, only minimal residual moisture should remain on the surface after the second wiping step. In the case of immersion disinfection, rinse the products off with water after the exposure time. Because of its good material compatibility, perform® can also be used on sensitive surfaces such as PMMA. Ensure that all visible dirt is removed before disinfection. Clean equipment thoroughly with water. For professional use only. Any serious incidents

related to the product are to be reported to the manufacturer and competent authority.

Information for order

Item	Delivery form	Item no.
perform® (60x) 40 g bag	60/Carton	on request
perform® (900 g) tin	4/Carton	on request
perform® (250x) 40 g bag	250/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Related Products

- mikrozid® PAA wipes
- terralin® PAA

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



GEPRÜFTES
UMWELTMANAGEMENT
REG. NR. 01-150-10003



ISO 9001:2008



ISO 13485

Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters

Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com