

schülke →



quartamon® med

Quaternary ammonium compounds-based concentrate for disinfecting and cleaning medical devices and surfaces.

Our Plus

- pleasant odour
- based on QAV

Application areas

quartamon® med is suitable for disinfecting surfaces and medical devices. As a long-established product, it is used in practice and hospital areas.

Instructions for use

Prepare the working solution in the desired concentration. Wet clean surfaces of medical devices, floors and other wipeable surfaces. Ensure that all visible dirt is removed before disinfection. Ensure complete wetting. Allows dosing via centralized and decentralized dosing devices. schülke dosing aids can be used for dosing. After completing the work, dispose of the working solution and clean the equipment thoroughly with water. For professional use only. Any serious incidents related to the product are to be reported to the manufacturer and competent authority. Not for final disinfection of semicritical and critical medical devices!

Microbiological efficacy

Efficacy	Concentration	Contact time
bactericidal <i>in accordance with VAH - clean conditions</i>	1 % (10 ml/l)	60 min.
bactericidal <i>in accordance with VAH - clean conditions</i>	2 % (20 ml/l)	15 min.
MRSA <i>EN13727 - clean conditions</i>	0.25 % (2.5 ml/l)	1 min.
yeastcidal <i>in accordance with VAH - clean conditions</i>	1 % (10 ml/l)	60 min.

Efficacy	Concentration	Contact time
yeastcidal <i>in accordance with VAH - clean conditions</i>	2 % (20 ml/l)	15 min.
virucidal against enveloped viruses <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline</i>	2 % (20 ml/l)	15 min.
Vaccinia virus <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - clean conditions</i>	0.5 % (5 ml/l)	15 min.
Vaccinia virus <i>in accordance with DVV (German Registered Association for Combating Viral Diseases)//RKI Guideline - dirty conditions</i>	2 % (20 ml/l)	15 min.
Polyoma SV40 <i>in accordance with DVV</i>	4 % (40 ml/l)	60 min.
Rotavirus <i>EN14476 - clean conditions</i>	1 % (10 ml/l)	5 min.

Certificates

- IHO listed
- VAH certificate

CE 0297



Product data

Composition:

100 g solution contains the following active substances: 9.5 g Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/ BKC (C12-C16))

Labelling according to Regulation (EC) No. 648/2004:

5 - 15 % nonionic surfactants, perfumes, benzyl salicylate, limonene, hexyl cinnamal

Chemical-physical data

Color	green
Density	ca. 1,01 g/cm ³ / 20 °C
Flash point	54 °C / Method : DIN 51755 Part 1
Form	liquid
pH	9,5 / 20 °C
Viscosity, dynamic	ca. 10 mPa*s / 20 °C / Method : DIN 54453

Special advice

Use disinfectants safely. Always read the label and product information before use.

Ensure that all visible dirt is removed before disinfection. Wearing protective gloves is advised. Keep away from heat and sources of ignition. For further information on safe handling and storage conditions, please refer to the safety data sheet. The shelf-life of the opened container can be found on the packaging.

Information for order

Item	Delivery form	Item no.
quartamon® med 2 l bottle	5/ Carton	on request
quartamon® med 5 l	1/ Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

Application aids

Application aids	Item no.
Can key for 5 + 10 l	135810
Decanter funnel	117901
measuring cup 500 ml	136101
measuring cup 50 ml	136102
Opener all-in-one (1l eurobottle/can/container)	70000804
schülke dosing feeder 5 l / 10 l (20ml / stroke)	117101
schülke tap for 5 l / 10 l can	135501

Related Products

- pursept® A Xpress

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: www.schuelke.com.

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: info@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com