tamaVet®

**Multi 2TA Test**

**(test cassette)**

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Rev. 2, 2019-12-17 RZ

### 1. Intended Use

The tamaVet® Multi 2TA Test is used for the detection of feline leukaemia virus antigen p27 (FeLV) and the detection of antibodies directed against the feline immunodeficiency virus (FIV) associated proteins p24 and gp41 in **whole blood, serum or plasma of cats**. The test is intended for use as an aid in the diagnosis of a single or mixed infection of FeLV and/or FIV and is designed for professional use only.

### 2. Introduction and Clinical Significance

The Feline Leukaemia Virus (FeLV) belongs to the family Retrovirus and the subfamily *Orthoretrovirinae*.

Infections with FeLV are viral diseases which occur frequently weaken the immune system of cats.

The feline immunodeficiency virus (FIV) belongs, like the human immunodeficiency virus (HIV), to the family of the retroviridae, genus lentivirus.

An FIV infection progresses through various stages similar to an HIV infection in humans: the acute stage, stable stage and terminal stage - also known as cat AIDS.

During the acute stage, which can last from several days to a several weeks, clinical signs are unspecific symptoms like fever, weight loss and diarrhoea. After the acute stage, the infected cat enters the stable stage, the cat seems to be clinically healthy and shows no symptoms. The stable stage can last for several years.

Symptoms within the final, terminal stage of infection vary depending on the different co-infections. The immune system of the cat is considerably weakened and prone to bacterial and viral infections. As there is no suitable therapy available up to now, FIV infections are almost always fatal.

### 3. Test Principle

The tamaVet® Multi 2TA Test is a rapid **sandwich immunoassay** developed for professional use by veterinarians. The assay employs specific antibodies/antigens which capture FeLV antigens or anti-FIV antibodies from the sample and thus visualise them. The tamaVet® Multi 2TA Test is a highly sensitive immunoassay and comprises two test strips (one for FeLV and another for FIV) housed in a handy test cassette.

The tamaVet® FeLV test strip contains a sample pad, a conjugate pad, a reaction membrane and an absorption pad. The conjugate pad contains specific gold-labeled anti-FeLV antibodies. The reaction membrane is precoated with additional anti-FeLV antibodies immobilised in the test line region and, as a control, antibodies immobilised in the control line region. The absorption pad soaks up the fluid applied to the test cassette.

Following the application of the sample onto the test strip, FeLV antigens in the sample material react with the gold-labeled antibodies to form a complex in the conjugate pad. The fluid then migrates along the strip by capillary action. In the test line region, the formed complex is captured by the immobilised capture antibodies, causing the formation of a red line. The presence of a red line in the test line region

indicates a positive result. If the sample does not contain FeLV antigens, no line will form in the test line region. Additionally, a red line must form in the control line region irrespective of the presence of FeLV antigens in the specimen. The control line serves as an internal control and indicates that sufficient sample fluid has been applied to the test strip and that proper membrane wicking has occurred.

The tamaVet® FIV test strip consists of a sample pad, a conjugate pad, a reaction membrane and an absorption pad. The conjugate pad contains gold-labeled FIV antigens. The reaction membrane is precoated with second FIV capture antigens immobilised in the test line region and, as a control, antibodies immobilised in the control line region. The absorption pad soaks up the fluid applied to the test cassette.

Following the application of the sample onto the test strip, anti-FIV antibodies in the sample material react with the gold-labeled antigens to form a complex in the conjugate pad. The fluid then migrates along the strip by capillary action. In the test line region, the formed complex is captured by the immobilised capture reagents, causing the formation of a red line. The presence of a red line in the test line region indicates a positive result. If the sample does not contain anti-FIV antibodies, no line will form in the test line region. Additionally, a red line must form in the control line region, irrespective of the presence of anti-FIV antibodies in the specimen. The control line serves as an internal control and indicates that sufficient sample fluid has been applied to the test strip and that proper membrane wicking has occurred.

### 4. Reagents and Materials Supplied (5 tests/kit)

- 5 tamaVet® Multi 2TA test cassettes
- 2.5 mL buffer solution
- 5 pipettes (single use)
- 1 package insert

### 5. Additional Materials Required

- Timer

### 6. Storage & Stability

The tamaVet® Multi 2TA Tests must be stored at room temperature (15-30°C). The test is stable until the expiry date printed on the foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

### 7. Warnings and Precautions

- For veterinary use only.
- For professional *in-vitro* diagnostic use only.
- For single use only.
- Store tests between 15-30°C. Do not freeze tests.
- Do not use tests beyond the expiry date.
- Do not use tests with damaged pouches.
- Use the test within 60 minutes of removing it from the pouch.
- Consider the test results as invalid read after more than 10 minutes.

- Use the stated amount of sample. An incorrect number of drops can lead to a false result.
- Place the test cassette on a smooth, flat surface while the test is being performed.
- Do not touch the reaction area.
- Use only the test components included in the kit. Do not mix or replace the test components with components of a different kit.
- Avoid cross contaminations by using a new sample tube and a new pipette for each sample.
- Sample materials are potentially infectious. During the test procedure, adhere to standard guidelines for the handling of potentially infectious materials and chemical reagents. The use of protective clothing (laboratory coat, gloves, eye protection) is recommended. Materials that have come into contact with the sample material should be disposed of in accordance with local regulations.

### 8. Specimen Collection and Preparation

The tamaVet® Multi 2TA Test was developed for the examination of **whole blood, serum or plasma of cats**.

#### General comments:


- Perform the test immediately after specimen collection. If the test is not to be performed immediately after sampling, serum and plasma specimens may be stored for up to 3 days and whole blood samples for up to 2 days in a closed container at 2-8°C. For long-term storage, serum and plasma samples must be kept below -20°C. The samples should not be frozen and thawed repeatedly. Do not freeze whole blood specimens.
- Use only clear, non-hemolyzed serum or plasma specimens. **Heparin or EDTA blood** are also suitable for testing.
- Bring samples to room temperature prior to testing, though do not leave them at this temperature for prolonged periods of time.
- Testing whole blood may decrease the sensitivity of the test. If the test result is negative despite an existing suspicion of an infection, the test should be repeated with serum or plasma samples.
- Shake whole blood samples well before testing.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiological agents.
- The transportation of whole blood samples should only be performed in exceptional cases. The shipping of serum and plasma is preferable.

### 9. Test Procedure

1. Specimen and kit components must be brought to room temperature before performing the test.
2. Open the pouch and place the test cassette on a horizontal surface. Remove the pipette from the pouch and take up the sample material using the pipette.


#### 3a. Serum, Plasma

Holding the pipette vertically, dispense **1 drop (30µL) of serum or plasma** into each sample well and wait a couple of seconds until everything is absorbed.


#### 3b. Whole blood

Holding the pipette vertically, dispense **1 drop (30µL) of whole blood** into each sample well and wait a moment until the fluid is absorbed.

4. Open the buffer tube and, holding it vertically, add **2 drops of the buffer solution** to each sample well. Avoid the formation of air bubbles. Start the timer. If the liquid is not migrating sufficiently along the strip after 60 seconds for serum and plasma, or 90 seconds for whole blood, add an additional drop of the buffer fluid into the respective sample well.

5. Read the results 10 minutes after starting the timer. **Do not interpret the result after more than 10 minutes.** The test and the result are considered as invalid after this period.

### 10. Result Interpretation

#### Positive result

Each test is positive if the respective control line (C) and the test line (T) are visible in the reaction area. FeLV antigens and/or anti-FIV antibodies have been detected in the sample material. If a weakly defined test line appears, the test result is still positive.

#### Negative result

Each test is negative if only the control line (C) is visible in the reaction area. No FeLV antigens or anti-FIV antibodies have been detected in the sample material.

#### Invalid result

If no control line is visible after the test has been performed, the test is invalid. All results that show no control line after the denoted time must be discontinued. Review the test procedure and repeat the test with a new test cassette.

If problems persist, discontinue using the test kit and contact your distributor.

**Possible test results:**

- Both tests are positive (C- and T-lines appear in both reaction areas)
- Both tests are negative (only the C-lines appear)
- One test is positive (besides both C-lines only one T-line in one reaction area appears)

**11. Quality Control**

Several internal controls are included in the test cassette:

A red line forming in the control line region (C) is considered as internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. In addition, a green/blue line is visible in the control line region before starting the test. It is considered a quality control method and will be washed away by the sample liquid during the test.

**12. Limitations**

- Test results must always be interpreted in light of all the medical findings and never individually. All test results which do not correlate with the clinical picture must be repeated. Positive results detected with the tamaVet® veterinary rapid tests must be confirmed using another analytical method.
- There is a possibility that the test result can be distorted by technical errors, errors in the test procedure or by substances or factors that influence the test and have not been mentioned here.

**13. Performance Characteristics**

	Sensitivity	Specificity	OA
FeLV	94.6%	99.9%	97.2%
FIV	91.7%	97.8%	99.1%

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nal von minden GmbH