

SURGICAL GOWN «STANDARD»



DESCRIPTION

The **Teqler "Standard" surgical gown** is made from liquid repellent, breathable non-woven material (SMMS). With help of Velcro strips in the neck area and straps, the gown can be adjusted individually. Available in different sizes.

PRODUCT DETAILS

- Standard surgical gown
- Individually sterile packed with label stickers on every package
- Material: SMMS (Spunbond - Meltblown - Meltblown - Spunbond)
- Ideal for operations with light to moderate liquid occurrence
- High resistance to germ transmigration and seepage
- Wide overlapping back section, stable tie bands
- Knitted arm cuffs, velcro closure in neck area
- Colour: blue
- Available in sizes M, L, XL and XXL
- 1 surgical gown, delivered with 2 sterile non-woven hand towels

Overview

Article Number:
T135280

Category:
Protective clothing OR



Teqler is a trademark of NetMed S.à.r.l.

📍 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

☎ +352 267149 09 📠 +352 267149 19 ✉ info@netmed.lu

Social Media

