

PETRI DISH



DESCRIPTION

This petri dish is made of temperature and chemical-resistant borosilicate glass, making it suitable for use in the laboratory. The glass dish and the lid can both be autoclaved and are available in various sizes.

PRODUCT DETAILS

- Petri dish (transparent)
- Made of 3.3 borosilicate glass
- Resistant to high temperatures and many chemicals
- Can be autoclaved
- Available in the following sizes: 60 mm, 75 mm, 100 mm
- Delivered quantity: 1 piece




Overview

Article Number:
T135839-T135841

Category:
Lab Equipment

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

