

WASH BOTTLES WITH WHITE CAP | 250 ML



DESCRIPTION

This wash bottle is semi-transparent with a white screw cap made of durable plastic. The use of flexible low density polyethylene (LDPE) allows one to squeeze the wash bottle and optimally dose the liquid.

PRODUCT DETAILS

- Material (bottle): thermoplast / LDPE (low density polyethylene)
- Colour (bottle): semi-transparent / milky
- Material (cap): plastic
- Colour (cap): white
- Handy and easy to squeeze
- Size: 250 ml
- Delivered quantity: 1 piece

Overview

Article Number:
T135860

Category:
Lab Equipment

EAN:
4260547369398



Teqler is a trademark of NetMed S.à.r.l.

📍 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

☎ +352 267149 09 📠 +352 267149 19 ✉ info@netmed.lu

Social Media

