

## ALUMINIUM WHEELCHAIR



### DESCRIPTION

The robust Teqler **aluminium wheelchair** can easily be folded together and is equipped with push handles and hand brakes.

### PRODUCT DETAILS

- Lightweight, space-saving and smooth rolling
- Upholstered seat and backrest made of robust, low-maintenance material; colour: black
- Footrests can be swivelled to the side and completely removed
- Side sections can be folded up to ease entry by the patient
- High-quality, robust PU front wheels
- Main wheels with treaded tyres
- Hand brakes on each side
- Reflectors for good visibility
- Frame colour: silver-grey
- Weight: 14.7 kg
- Maximum load capacity: 120 kg

### DIMENSIONS

- Seat width: 46 cm
- Seat width: 41,5 cm
- Backrest height: 42 cm
- Handle height: 90 cm
- Diameter of main wheels: 59.5 cm
- Diameter of front wheels: 19 cm
- Folded dimensions: 90 x 78 x 28 cm (H x L x W)
- Calf length: 40 cm
- Front seat height: 48 cm
- Back seat height: 44 cm

### Overview

Article Number:  
**T135300**




Category:  
**Mobility and transport aids**

EAN:  
**4260306776474**



**Teqler is a trademark of NetMed S.à.r.l.**

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

**Social Media**

