



TEQLER MEDICAL

# MEDICAL

## TEQLER DISPOSABLE INJECTION NEEDLES



### DESCRIPTION

The **Teqler disposable injection cannulas** made of chrome-nickel steel have a colour-coded cannula attachment. The silicone coating ensures a painless injection. Available in different sizes.

### PRODUCT DETAILS

- With silicone coating for a painless injection
- Disposable cannula made of chrome-nickel steel
- Latex-, PVC- and DEHP-free
- With colour-coded plastic attachment
- Luer connector
- Complies with ISO 7864
- Individually sterile packed
- 1 package of 100 pieces

### Overview

Article Number:  
**T138333-T138340**

Category:  
**Infusion | Injection | Syringes**



Teqler is a trademark of NetMed S.à.r.l.

NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

+352 267149 09 +352 267149 19 info@netmed.lu

Social Media

