

EMR-validated flat scale with ID-Display and stable glass platform

seca
Precision for health

 Order here!

 **PRAXISDIENST**
Medical Supplies since 1953

seca 651

- + Display shows name or identification number (ID) of users and patients for enhanced safety
- + Data transmission controlled on the touchscreen
- + Large platform made of tested seca Bearclaw Glass® for absolute breakage resistance and perfect hygiene
- + Glass surface is easy to clean and disinfect (more hygienic than plastic)
- + EMR-validated: transmit measured data directly to any Electronic Medical Record (EMR) system

EMR



seca 651

EMR-validated flat scale with ID-Display and stable glass platform



Large, high-capacity flat scale developed for EMR integration

With the seca Scale-up Line, there is now the first flat scale developed from the start for integration with electronic medical records. It is not only EMR-validated, but also has an extremely large, low-profile platform which offers ample space even for particularly obese patients. The seca 650 from the seca Scale-up Line makes medical weighing for patients and users safer, more efficient and more comfortable.



seca Bearclaw Glass® platform

The glass platform not only looks good, it is also particularly hygienic. The death rate of many types of germs on glass is over 99 %. This makes the glass platform similarly hygienic as stainless steel and also easy to clean. Anti-slip print ensures that the glass is slip-resistant. Break-resistance is guaranteed by use of special seca Bearclaw Glass®.



ID-Display ensures an easy and user-friendly data transmission to EMR

seca Scale-up Line is developed especially for perfect integration of measured data with EMR systems. Weight and data such as time of measurement are stored centrally and are accessible for every user. The ID-Display shows name or patient ID with every measurement, which are gathered by an additional barcode scanner. This makes the integration process transparent, safe and even more user-friendly. For more information please visit integration.seca.com

| | |
|---------------------------------|---|
| Approval class | III |
| Capacity | 300 kg |
| Graduation | 50 g < 150 kg > 100 g |
| Dimensions platform (W x H x D) | 600 x 70 x 467 mm |
| Functions | Display of user and patient name or ID, customizable shortcut buttons, TARE, Pre-TARE, mother+child function, HOLD, Auto-HOLD, Auto-BMI, CLEAR, Auto-CLEAR, kg/lbs switch-over, cm/inch switch-over, damping, Wi-Fi, Ethernet |
| Accessories | Ultrasonic height measurement seca 257 in combination with adapter seca 464, Handrail seca 455 |

We reserve the right to make modifications.

seca united kingdom
 40 Barn Street / Digbeth T +44 121 643 9349
 Birmingham B5 5QB F +44 121 633 3403
 United Kingdom E info.uk@seca.com

Precision for health
seca.com/scaleup

20111102420