

TOP-NOTCH
PERFORMER

iFuge M15K
(Table Top Centrifuge)

iFuge M15K

(Table Top Centrifuge)

In any research facility or lab where accuracy reigns supreme, you need an equipment that can deliver precision to the last decimal point. Such an equipment might be hard to get hold of but not impossible, especially when Neutation is here to help you.

With Neutation's iFuge M15K you can expect a lot as it delivers even more. The lightning fast 15000 RPM microprocessor controlled centrifuge can perform at peak efficiency delivering top-notch accuracy in every cycle despite high load.




Witness The Champion



Brushless Dc Motor For
Maintenance Free Long Life



Variable Speed
From 500-15000 RPM



Microprocessor
Controlled Efficient
Performance.



Intuitive Interface For
Quick And Convenient
Setting With Last Run
Memory Feature



One Touch RPM / RCF
Conversion Operation

Variable Timer Setting
From 1 To 99 Min &
Infinite Mode



One-touch
Short Spin Feature



Easy fit rotor



12 X 1.5 / 2 ml
Closed Rotors Design
For Lesser Heat Generation



Imbalance Detection
With Auto-cut Off
Safety Feature



Witness The Champion



Lid-lock Safety Feature
Opens The Lid Automatically
After Every Run



Enhanced Airflow Design
For Less Noise And
Greater Sample Safety



Smooth & Silent Operation
Even At Maximum Speed

Safe Low
Voltage Device



Compact Design And
Can Also Be Used In
Hoods And Cold Rooms



ROTOR OPTION



12 x 1.5/2 ml | 15596 x g



2 x 8 x 0.2 ml strip rotor

Specifications

iFuge M15K	
Specifications	Details
Rotor Capacity	12 x 1.5ml / 2.0ml Microtubes
Maximum RPM / RCF	15000 / 15596g
Display	Large Digital
Imbalance Detection	Yes
Motor Type	(BLDC) Brushless DC motor
Timer Setting	Countdown Timer From 1-99 Mins & Infinity Mode
Dimensions (WxDxH)	187 x 265 x 115 mm
Noise Level	< 60 db
Weight	2.2 Kg
Rotors and Adaptors (included in standard pack)	Reduction Adaptors For 0.2 ml Microtubes, Reduction Adaptors For 0.5 ml Microtubes

Ordering Information

iFuge M15K		
Product Name	Cat #	Details
iFuge M15K	M112-5115 UK/AU/US/EU	High Speed Micro Centrifuge Upto 15000 RPM

Accessories		
Part Name	Cat #	Details
12 Slot Rotor	M12-R12	12 x 1.5/ 2 ml tubes slot rotor
Rotor Lid	M12-R12E	Rotor Lid for 12 slot rotor
Rotor Nut	M12-RN	Rotor Nut for the Rotor Lid
PCR Strip Rotor	M12 - R2 x 8	2 PCR stirp rotor in Black colour
PCR Strip Rotor Lid	M12-R2 x 8E	-
Adaptors for 0.2 ml	M12-D1	Reduction Adaptors of 0.2 ml tubes
Adaptors for 0.4/0.5 ml	M12-D2	Reduction Adaptors of 0.4/0.5 ml tubes
Power Supply adaptor	M12-A12	Power adaptor with specified pin type

NEUATION
TECHNOLOGIES
www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

*We also execute OEM projects. Specifications can be changed without notice for quality improvement.