


# SUPRAMID

BLACK



<b>Type</b>	Twisted and coated multifilament Sizes 5/0 and 6/0 monofilament
<b>Composition</b>	Polyamide 6.6 strands enclosed in a polyamide 6 sheet
<b>Colour</b>	Black
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	Non absorbable, gradually encapsulated by connective tissue
<b>Presentation</b>	Box with 12 sutures; Cassette

**PRAXISDIENST**  
Medical Supplies since 1953  
**Order here!**







*For quality  
&  
safety*



**smi**



[www.sutures.be](http://www.sutures.be)






 Round bodied taper point		Metric	USP	Length	Order number
<b>1/2 circle</b> 	HR 16 mm	1.5	4/0	75 cm	<i>black</i> 4150116
		2	3/0	75 cm	<i>black</i> 4200116
	HR 20 mm	2	3/0	45 cm	<i>black</i> 4200120
		3	2/0	75 cm	<i>black</i> 4300126
	HR 26 mm	3	2/0	75 cm	<i>black</i> 4300126
		3	2/0	75 cm	<i>black</i> 4300126


 Round bodied taper point		Metric	USP	Length	Order number
<b>Straight</b> 	GR 22 mm	1.5	4/0	50 cm	<i>black</i> 4152022


 Reverse cutting		Metric	USP	Length	Order number
<b>1/2 circle</b> 	HS 17 mm	1.5	4/0	75 cm	<i>black</i> 4151417
		2	3/0	75 cm	<i>black</i> 4201417
	HS 20 mm	1	5/0	75 cm	<i>black</i> 4101420
		1.5	4/0	75 cm	<i>black</i> 4151420
	HS 22 mm	2	3/0	75 cm	<i>black</i> 4201420
		1.5	4/0	75 cm	<i>black</i> 4151422
	HS 26 mm	2	3/0	75 cm	<i>black</i> 4201426
		3	2/0	75 cm	<i>black</i> 4301426

 Cutting		Metric	USP	Length	Order number
<b>Straight</b> 	GS 20 mm	1.5	4/0	50 cm	<i>black</i> 4152120
	GS 25 mm	1.5	4/0	60 cm	<i>black</i> 4152125

 Reverse cutting		Metric	USP	Length	Order number
<b>3/8 circle</b> 	DS 12 mm	0.7	6/0	75 cm	<i>black</i> 4071512
		1	5/0	45 cm	<i>black</i> 4101512-45
		1	5/0	75 cm	<i>black</i> 4101512
		1.5	4/0	45 cm	<i>black</i> 4151512

Reverse cutting		Metric	USP	Length	Order number
 3/8 circle	DS 16 mm	0.7	6/0	75 cm	<i>black</i> 4071516
		1	5/0	45 cm	<i>black</i> 4101516-45
		1	5/0	75 cm	<i>black</i> 4101516
		1.5	4/0	45 cm	<i>black</i> 4151516
		1.5	4/0	75 cm	<i>black</i> 4151516-75
		2	3/0	45 cm	<i>black</i> 4201516
		2	3/0	75 cm	<i>black</i> 4201516-75
	DS 16 mm black	1	5/0	75 cm	<i>black</i> 4101516VB
	DS 19 mm	0.7	6/0	75 cm	<i>black</i> 4071519
		1	5/0	75 cm	<i>black</i> 4101519
		1.5	4/0	45 cm	<i>black</i> 4151519
		1.5	4/0	75 cm	<i>black</i> 4151519-75
		2	3/0	45 cm	<i>black</i> 4201519-45
		2	3/0	75 cm	<i>black</i> 4201519
DS 19 mm Easy Pass	1.5	4/0	75 cm	<i>black</i> 4151519P-75	
DS 19 mm black	1	5/0	75 cm	<i>black</i> 4101519VB	
	1.5	4/0	75 cm	<i>black</i> 4151519VB	
DS 24 mm	2	3/0	75 cm	<i>black</i> 4201524	
	3	2/0	75 cm	<i>black</i> 4301524	
	3.5	0	75 cm	<i>black</i> 4351524	
DS 36 mm	3	2/0	75 cm	<i>black</i> 4301536	
	3.5	0	75 cm	<i>black</i> 4351536	
	4	1	75 cm	<i>black</i> 4401536	

Round bodied taper cutting point		Metric	USP	Length	Order number
 1/2 circle	HRT 18 mm	1.5	4/0	75 cm	<i>black</i> 4150518
		2	3/0	75 cm	<i>black</i> 4200518

Round bodied taper cutting point		Metric	USP	Length	Order number
 3/8 circle	DRT 18 mm	1.5	4/0	75 cm	<i>black</i> 4150418
		2	3/0	75 cm	<i>black</i> 4200418

### IN CASSETTES

Metric	USP	Length	Order number
1.5	4/0	100 m	4015
2	3/0	100 m	4020
3	2/0	100 m	4030
3.5	0	100 m	4035
4	1	100 m	4040
5	2	75 m	4050
6	3&4	50 m	4060
7	5	40 m	4070
8	6	40 m	4080
9	7	30 m	4090



**PRAXISDIENST**  
Medical Supplies since 1953

**Order here!**