

# Veterinary Automatic Air Warming System

## 01 Product Introduction


WARMAIR for Veterinary Automatic Air Warming System (hereinafter referred to as WARMAIR) consists of a five-range temperature management module, an internal cleaning and disinfection module, a hot air purification module.

WARMAIR is connected to the warming air cushion through a hose. Warm air is generated in the main unit and forced into the air cushion through the hose. The warming air cushion of a proper model is placed around, on or under the pet to receive a surgery. Warm air will be evenly blown onto the pet through the small holes on the warming air cushion.




## 02 Scope of Application

WARMAIR for Veterinary Automatic Air Warming System is mainly used to prevent and treat hypothermia. In addition, if the environment conditions may cause drop in the body temperature of the sick pet, WARMAIR can be used to provide a comfortable temperature for the sick pet. WARMAIR is suitable for both adult pets and young pets.



**Smart  
Temperature Control**



**TOOTOO MEDITECH CO., LTD.**

Address: Office Room 607, Huate Commercial Building, NanLian Rd, Longgang District, ShenZhen City, China.

Tel: 86(755)2830 6385


Web: [www.newtootoo.com](http://www.newtootoo.com)



# Warm Air

VETERINARY AUTOMATIC AIR WARMING SYSTEM

## Veterinary Automatic Air Warming System



Air cushion



Hospital cage adaptor



Hose + Leather sheath

**TOOTOO MEDITECH CO., LTD.**

Address: Office Room 607, Huate Commercial Building, NanLian Rd, Longgang District, ShenZhen City, China.

Tel: 86(755)2830 6385

Web: [www.newtootoo.com](http://www.newtootoo.com)