

INJECTION PLASTER



DESCRIPTION

The **Teqler injection plasters** are available in two sizes and suitable even for sensitive skin.

PRODUCT DETAILS

- Wound cover made of non-woven material
- Does not adhere to the wound
- Also suitable for sensitive skin
- Non-sterile
- Skin coloured
- Available in sizes 4 x 1 cm and 6 x 2 cm
- 1 box with 400 pieces (4 x 1 cm) or with 200 pieces (6 x 2 cm)

Overview




Article Number:
T603930-T603931

Category:
**Infusion | Injection |
Syringes**

Packaging:
**1 box of 400 pieces (4 x 1
cm) or 200 pieces (6 x 2
cm)**

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

