



SUZHOU SUNING UNDERPAD CO., LTD
蘇州市蘇寧床墊有限公司

Technical data sheet

Item No.:	820501
Product name:	Baby underpad
Product description:	Baby underpad for incontinence usage, mostly in hospital, homecare ,babycare and so on.
Sizes of pad:	40 x 60 cm
Inner absorbent pad size:	40 x 50 cm
Raw Materials:	Nonwoven,Tissue, PE,Fluff pulp,SAP
Total weight per piece:	40 g (14.35 g Amount of top nonwoven+tissue+ PE) (23.65 g Amount of fluff pulp) (2 g Amount of SAP) (22 gms Back layer PE) (13 gms Tissue) (12 gms Top layer nonwoven)
Absorption:	0.5 liters absorption capacity
Color:	Purple
Packaging:	10 pcs / bag 5 bags / carton
Packing Materials:	Polyethylene 40 µm thickness
Carton size:	29 x 27 x 22 cm
Useful Life:	5 years
Handling and Storage:	Store dry,Avoid direct sunlight
Storage temperature:	temperature between -5-40 °C

Cautions, Warnings or contraindications

No. 70 Yong'an Road, Xuguan Industrial Park, High-Tech District, Suzhou, Jiangsu, China 215151
Telephone: 86-512-66653331, Fax: 86-512-66651066
www.sn-anpa.com





SUZHOU SUNING UNDERPAD CO., LTD
蘇州市蘇寧床墊有限公司

Not to use if the product is damaged or dirty, or do have other functional failures.

Rubbish and Final Disposition

The product must be disposed, disposable waste regulation, after using.

Date: 2021-11-12