

## OXYGEN NASAL CANNULA



### DESCRIPTION

The Teqler nasal cannula is DEHP and DOP-free and equipped with a standard connector. Nasal cannulas are used for nasal oxygen administration, e.g. in symptomatic therapy of an COPD exacerbation. While oxygen therapy can be performed with an oxygen mask or nasal cannula, the latter is preferable as it offers a higher level of patient comfort.

The tube brackets are positioned at the ears and attached at the patients chin with help of straps to prevent the cannula from moving. The prongs are placed in the nostrils. The Teqler nasal cannula is delivered with a 2.10m long tube.

### PRODUCT DETAILS

- Nasal cannula for adults
- For short and long-term oxygen therapy (e.g. in COPS or pulmonary emphysema)
- DEHP and DOP-free
- Latex-free
- Nasal cannula with standard connector
- Non-sterile
- Tube length: 2.10m
- PZN: 3048929
- 1pc.

### Overview

Article Number:  
**T999229**

Category:  
**General Supplies**



EAN:  
**4034342010218**

Dimensions:  
**Schlauch-Länge der  
Sauerstoff-Brille: 2,10 m**



**Teqler is a trademark of NetMed S.à.r.l.**

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

**Social Media**

