

**AeroDawg** 

Canine Aerosol Chamber

SMALL

MEDIUM

LARGE

## Mask Sizer

Each device is shipped with  
2 different sized masks.

**AeroDawg\*** Chamber SMALL,  
includes small and medium masks

**AeroDawg\*** Chamber LARGE,  
includes medium and large masks

## How to Use

The mask selected should cover the nose and mouth. However, if not possible  
apply mask to the dog's nostrils while gently holding the muzzle closed.



Trudell Medical International\*

[www.aerodawg.com](http://www.aerodawg.com)



201092-001 Rev.A. \* trade-marks and registered trade-marks of Trudell Medical International.  
© TMI 2015. All rights reserved.